

BetteStarlet



Information




The BetteStarlet, 1900 × 900 × 420 mm, article no. 1830, is not available with BetteUpstand.

For professional installation we recommend the use of the original bath feet B23-1500.

If the bath is to be undermounted, template information can be provided.

Scope of delivery

Including noise-reducing pads.



Dimensions

article no.	dimensions in mm	A	B	C	D	E	F	G	H	I	K	T	U	V	volume*
1210	1500 × 800 × 420	1500	800	80	1340	80	540	420	540	1070	555	495	80	80	107 litre
2590	1570 × 700 × 420	1570	700	65	1440	65	470	630	470	1170	455	425	65	95	92 litre
2540	1600 × 650 × 420	1600	650	80	1440	80	480	640	480	1170	455	425	55	55	92 litre
1220	1600 × 700 × 420	1600	700	80	1440	80	485	630	485	1170	455	425	80	80	92 litre
1230	1650 × 700 × 420	1650	700	105	1440	105	510	630	510	1170	455	425	80	80	92 litre
1730	1700 × 700 × 420	1700	700	80	1540	80	505	690	505	1270	505	480	55	55	122 litre
1380	1700 × 750 × 420	1700	750	80	1540	80	500	700	500	1270	505	480	80	80	122 litre
1450	1750 × 800 × 420	1750	800	105	1540	105	520	710	520	1270	505	485	105	105	122 litre
1430	1800 × 750 × 420	1800	750	80	1640	80	530	740	530	1370	555	485	55	55	155 litre
1630	1800 × 800 × 420	1800	800	80	1640	80	530	740	530	1370	555	485	80	80	155 litre
1631	1850 × 850 × 420	1850	850	105	1640	105	562	726	562	1370	555	485	105	105	155 litre
1830	1900 × 900 × 420	1900	900	85	1730	85	560	780	560	1465	625	560	80	80	205 litre

*Water content less 70 litres (for human displacement).